



## Analytical and Clinical Validation of Expressed Variants and Fusions From the Whole Transcriptome of Thyroid FNA Samples

Trevor E. Angell<sup>1\*</sup>, Lori J. Wirth<sup>2</sup>, Maria E. Cabanillas<sup>3</sup>, Maisie L. Shindo<sup>4</sup>, Edmund S. Cibas<sup>5</sup>, Joshua E. Babiarz<sup>6</sup>, Yangyang Hao<sup>6</sup>, Su Yeon Kim<sup>6</sup>, P. Sean Walsh<sup>6</sup>, Jing Huang<sup>6</sup>, Richard T. Kloos<sup>7</sup>, Giulia C. Kennedy<sup>8</sup> and Steven G. Waguespack<sup>3</sup>

<sup>1</sup> Division of Endocrinology, Diabetes and Metabolism, Keck School of Medicine, University of Southern California, Los Angeles, CA, United States, <sup>2</sup> Department of Medicine, Massachusetts General Hospital, Boston, MA, United States, <sup>3</sup> Department of Endocrine Neoplasia and Hormonal Disorders, The University of Texas MD Anderson Cancer Center, Houston, TX, United States, <sup>4</sup> Otolaryngology–Head & Neck Surgery, Oregon Health & Science University, Portland, OR, United States, <sup>5</sup> Department of Pathology, Brigham and Women's Hospital and Harvard Medical School, Boston, MA, United States, <sup>6</sup> Research and Development, Veracyte, South San Francisco, CA, United States, <sup>7</sup> Medical Affairs, Veracyte, South San Francisco, CA, United States, <sup>8</sup> Research and Development, Medical Affairs, and Clinical Affairs, Veracyte, South San Francisco, CA, United States

#### **OPEN ACCESS**

#### Edited by:

Vasyl Vasko, Uniformed Services University of the Health Sciences, United States

#### Reviewed by:

Joanna Klubo-Gwiezdzinska, National Institutes of Health (NIH), United States Dorina Ylli, MedStar Health Research Institute (MHRI), United States

> \*Correspondence: Trevor E. Angell trevor.angell@med.usc.edu

#### Specialty section:

This article was submitted to Thyroid Endocrinology, a section of the journal Frontiers in Endocrinology

Received: 01 May 2019 Accepted: 22 August 2019 Published: 11 September 2019

#### Citation:

Angell TE, Wirth LJ, Cabanillas ME, Shindo ML, Cibas ES, Babiarz JE, Hao Y, Kim SY, Walsh PS, Huang J, Kloos RT, Kennedy GC and Waguespack SG (2019) Analytical and Clinical Validation of Expressed Variants and Fusions From the Whole Transcriptome of Thyroid FNA Samples. Front. Endocrinol. 10:612. doi: 10.3389/fendo.2019.00612 **Introduction:** The Afirma<sup>®</sup> Xpression Atlas (XA) detects gene variants and fusions in thyroid nodule FNA samples from a curated panel of 511 genes using whole-transcriptome RNA-sequencing. Its intended use is among cytologically indeterminate nodules that are Afirma GSC suspicious, Bethesda V/VI nodules, or known thyroid metastases. Here we report its analytical and clinical validation.

**Methods:** DNA and RNA were purified from the same sample across 943 blinded FNAs and compared by multiple methodologies, including whole-transcriptome RNA-seq, targeted RNA-seq, and targeted DNA-seq. An additional 695 blinded FNAs were used to define performance for fusions between whole-transcriptome RNA-seq and targeted RNA-seq. We quantified the reproducibility of the whole-transcriptome RNA-seq assay across laboratories and reagent lots. Finally, variants and fusions were compared to histopathology results.

**Results:** Of variants detected in DNA at 5 or 20% variant allele frequency, 74 and 88% were also detected by XA, respectively. XA variant detection was 89% when compared to an alternative RNA-based detection method. Low levels of expression of the DNA allele carrying the variant, compared with the wild-type allele, was found in some variants not detected by XA. 82% of gene fusions detected in a targeted RNA fusion assay were detected by XA. Conversely, nearly all variants or fusions detected by XA were confirmed by an alternative method. Analytical validation studies demonstrated high intra-plate reproducibility (89%-94%), inter-plate reproducibility (86–91%), and inter-lab accuracy (90%). Multiple variants and fusions previously described across the spectrum of thyroid cancers were identified by XA, including some with approved or investigational targeted therapies. Among 190 Bethesda III/IV nodules, the sensitivity of XA as a standalone test was 49%.

1

Conclusion: When the Afirma Genomic Sequencing Classifier (GSC) is used first among Bethesda III/IV nodules as a rule-out test, XA supplements genomic insight among those that are GSC suspicious. Our data clinically and analytically validate XA for use among GSC suspicious, or Bethesda V/VI nodules. Genomic information provided by XA may inform clinical decision-making with precision medicine insights across a broad range of FNA sample types encountered in the care of patients with thyroid nodules and thyroid cancer.

Keywords: molecular diagnostics, thyroid cancer, fine-needle aspiration, thyroid molecular assays, RNAsequencing, transcriptome, atypia of undetermined significance, follicular neoplasm

#### INTRODUCTION

Genomic assessment for precision medicine is a story of incredible advancement. In the Nineteenth century, observations of dividing cancer cells suggested that they were abnormal clones caused by defects of hereditary material (1). The first isolation of a specific DNA variant responsible for cancer formation was in 1982, a G>T substitution in codon 12 of the HRAS gene (1). Since then there has been an explosion of cancer genome understanding, with the documentation of more than 350 cancer driving genes and 100,000 somatic mutations by the early Twenty-first century (1). Subsequently, a common theme that has emerged in oncology suggests that each cancer can be genomically subtyped and that the downstream gene expression profile predicts its cellular morphology, clinical presentation, prognosis, and that this presents an opportunity for the development of effective targeted therapies.

Investigation and understanding of benign and malignant thyroid nodules has followed a similar course as many insights have been gained from large genomic studies across a spectrum of histologic subtypes (2-14). Concomitantly, targeted therapies for advanced thyroid cancer have emerged, including FDA approved or investigational selective inhibitors of AKT<sup>1</sup>, ALK (15, 16), BRAF<sup>1</sup> (16, 17), cKIT<sup>1</sup>, EGFR<sup>2</sup>, HRAS<sup>3,4,5,6</sup> (18), KRAS<sup>2,4,5,6,7,8</sup>, MET<sup>1</sup>, mTOR (16, 19), NRAS<sup>1,4,5,6</sup>, NTRK (15, 16, 20), PAX8/PPARG<sup>9</sup>, PIK3CA (PI3K) (21), PTEN<sup>1</sup>, RET (16, 22, 23), ROS1 (15, 16), and microsatellite instability-high or mismatch repair deficient solid tumors (24).

Currently, an important tool to avoid unnecessary diagnostic surgery among cytologically indeterminate thyroid nodules is the Afirma Genomic Sequencing Classifier (GSC) (25, 26). In this context, "cytologically indeterminate" refers to the two Bethesda categories Atypia of Undetermined Significance/ Follicular Lesion of Undetermined Significance ("Bethesda III") and Follicular Neoplasm/ Suspicious for a Follicular Neoplasm

(Bethesda IV) (27) or their equivalents (28). The Afirma GSC is a cancer rule-out test that partners whole transcriptome RNA sequencing genomic information derived from a fineneedle aspiration (FNA) biopsy with machine learning to create algorithms that identify specific neoplasms, including MTC (see Table 1 for histology subtype abbreviations), and ultimately classify the sample as GSC benign or suspicious. Nodules identified as GSC benign have a cancer risk of approximately 4% and can be considered for clinical observation in lieu of diagnostic surgery (29-31). Conversely, GSC suspicious nodules have an increased cancer risk of approximately 50%, which is roughly 2-fold higher than it was based on cytology alone. These nodules are typically considered for surgical resection (29-31). The GSC algorithms rely heavily on differential gene expression for sample classification. Several of the included modules make limited use of RNA-sequencing's ability to detect genomic variants and fusions in the transcribed RNA, including BRAF V600E variants, and RET/PTC1 and RET/PTC3 fusions, which are all highly predictive of malignancy. However, it was not until the introduction of the Afirma Xpression Atlas (XA) that variant and fusion identification by RNA-sequencing was more significantly harnessed. The use of whole transcriptome sequencing by both the Afirma GSC and Afirma XA allows the same sample collection and shipping method to be used for both tests, and both tests are run on the same FNA sample. Using one sample for both tests facilitates successful test results despite the small genomic sample obtained by FNA. XA findings may predict tissue cellular morphology, clinical syndromes, cancer behavior (including mode of metastasis), prognosis, and facilitate the selection of effective targeted therapy in the appropriate clinical setting.

Here we describe the analytical and clinical validation of XA to report nucleotide variants and gene fusions beyond BRAF V600E, RET/PTC1, and RET/PTC3 fusions using whole transcriptome RNA-seq data derived from FNA samples. We compared the XA results to a targeted DNA panel for nucleotide variants and compared XA results to a targeted RNA fusion panel for gene/gene fusion detection. The data demonstrate a high level of agreement between methods and that these variant and fusion methods alone cannot serve as rule-out test to exclude cancer/noninvasive follicular thyroid neoplasm with papillary-like nuclear features (NIFTP). When used among Bethesda III/IV nodules that are GSC suspicious, Bethesda V/VI thyroid nodules, or by extension known thyroid cancer

<sup>&</sup>lt;sup>1</sup>https://clinicaltrials.gov/ct2/show/NCT02465060

<sup>&</sup>lt;sup>2</sup>https://clinicaltrials.gov/ct2/show/NCT03065387

<sup>&</sup>lt;sup>3</sup>https://clinicaltrials.gov/ct2/show/NCT02383927 <sup>4</sup>https://clinicaltrials.gov/ct2/show/NCT03244956

<sup>&</sup>lt;sup>5</sup>https://clinicaltrials.gov/ct2/show/NCT03181100 <sup>6</sup>https://clinicaltrials.gov/ct2/show/NCT00019331

<sup>&</sup>lt;sup>7</sup>https://clinicaltrials.gov/ct2/show/NCT03600883

<sup>&</sup>lt;sup>8</sup>https://clinicaltrials.gov/ct2/show/NCT03785249

<sup>9</sup>https://clinicaltrials.gov/ct2/show/NCT01655719

#### TABLE 1 | Histopathology subtypes.

#### Malignant subtypes

malignant su	btypes
FC	Follicular carcinoma. Variants include capsular invasion (FC-c) and vascular invasion (FC-v) $$
FVPTC	Follicular variant of papillary thyroid carcinoma. Variants include FVPTC micro carcinomas (mFVPTC)
HCC	Hürthle cell carcinoma. Variants include capsular invasion (HCC-c) and vascular invasion (HCC-v) $$
PTC	Papillary thyroid carcinoma. Variants include PTC micro carcinomas (mPTC), tall-cell variant (PTC-TCV), and tall-cell variant micro carcinomas (mPTC-TCV)
MTC	Medullary thyroid carcinoma
PDC	Poorly differentiated carcinoma
WDC-NOS	Well-differentiated carcinoma not otherwise specified
Benign subty	
Benign subty	pes
Benign subty BFN	<b>pes</b> Benign follicular nodule Chronic lymphocytic thyroiditis (aka, Hashimoto's thyroiditis). Also
Benign subty BFN CLT	<b>pes</b> Benign follicular nodule Chronic lymphocytic thyroiditis (aka, Hashimoto's thyroiditis). Also known as LCT (lymphocytic thyroiditis)
Benign subty BFN CLT FA	<b>pes</b> Benign follicular nodule Chronic lymphocytic thyroiditis (aka, Hashimoto's thyroiditis). Also known as LCT (lymphocytic thyroiditis) Follicular adenoma
Benign subty BFN CLT FA HCA	pes Benign follicular nodule Chronic lymphocytic thyroiditis (aka, Hashimoto's thyroiditis). Also known as LCT (lymphocytic thyroiditis) Follicular adenoma Hürthle cell adenoma
Benign subty BFN CLT FA HCA HN	pes Benign follicular nodule Chronic lymphocytic thyroiditis (aka, Hashimoto's thyroiditis). Also known as LCT (lymphocytic thyroiditis) Follicular adenoma Hürthle cell adenoma Hyperplastic nodule

metastases (data not shown), the genomic content provided by XA can provide additional information that may inform clinical decision-making.

#### MATERIALS AND METHODS

#### **Nucleic Acid Extraction**

RNA and DNA were extracted from thyroid FNAs or control tissues (described below under Cohorts and Controls, respectively) using the Qiagen AllPrep Micro Kit (Qiagen, Hilden, Germany) according to manufacturer's instructions. RNA was quantitated using Quantiflour (Promega, Madison, WI) and DNA was quantitated using PicoGreen (Promega, Madison, WI). Fluorescence was read on a Tecan Infinite M200 Pro (Tecan, Männedorf, Switzerland). RNA Integrity Number (RIN) was determined for RNA using RNA Pico Chips on the Bioanalyzer 2100 (Agilent, Santa Clara, CA).

#### **Custom AmpliSeq Panels**

The Ion AmpliSeq Designer (Thermo Fisher, Waltham, MA) was used to generate a custom DNA panel against 568 targets described in Pagan et al. (6). Nine Y chromosome SNPs were included to assign genomic gender, which was compared to clinical gender to ensure sample identity. A custom RNA AmpliSeq Variant Panel was designed against the same targets described for the DNA AmpliSeq panel. A custom AmpliSeq RNA Fusion Panel was also generated (Thermo Fisher, Waltham, MA), with 168 fusions plus 6 house-keeping genes for controls.

#### AmpliSeq Library Preparation

Ion Torrent Libraries were generated using the Ion AmpliSeq Library Kit 2.0 (Thermo Fisher, Waltham, MA) according

to manufacturer's instructions using 10 ng input material. RNA was first reverse transcribed with SuperScript IV VILO (Thermo Fisher, Waltham, MA). Following cleanup, libraries were quantitated using qPCR on a QuantStudio 6 (Thermo Fisher, Waltham, MA), normalized to 100 pM, and loaded onto the IonChef (Thermo Fisher, Waltham, MA). Pooled libraries were loaded onto the Ion 540 Chip (Thermo Fisher, Waltham, MA) and sequenced on the Ion Torrent S5XL (Thermo Fisher, Waltham, MA).

#### Controls

Each DNA plate included the Horizon Quantitative Multiplex Reference Standard (HD701, Horizon, Cambridge, United Kingdom), NA12878 (Coriell, Camden, NJ), *BRAF* V600E positive thyroid tissue (Cooperative Human Tissue Network, National Cancer Institute, Bethesda, MD), and *TERT* C228T positive thyroid tissue (Asterand, Westbury, NY).

## Targeted DNA/RNA Sequencing Data Analysis

The Ion Torrent Suite Software version 5.6.0 (Thermo Fisher, Waltham, MA) was used to demultiplex, map reads to the reference genome hg19, and call variants. Specifically, tmap version 5.6.8 was used for reads mapping, and Torrent Variant Caller version 5.6–10 was used to detect variants with parameter settings optimized for low frequency variant detection with minimal false negative calls on Ion AmpliSeq experiments.

### Targeted RNA Sequencing for Fusion Data Analysis

The Ion Reporter version 5.6.0 was used to detect fusions. Key fusion detection parameters were set as follows: minimum 20,000 total valid mapped reads to qualify a sample for further analysis; minimum 20 reads required to call a fusion; medium sensitivity, which requires 70% overlap between reads and reference sequence with at-least 66.66% exact matches in the overlap.

#### qPCR Data

Taqman assays were obtained from Thermo Fisher or designed and synthesized by Integrated DNA Technologies (IDT, Coralville, IA). Ten nanogram of RNA was reverse transcribed with QuantiTect (Qiagen, Hilden, Germany). qPCR reactions were performed in duplicate using Fast Advanced Master Mix (Thermo Fisher, Waltham, MA), VIC-labeled primer-limited TBP (Hs00427620\_m1, Thermo Fisher), and the FAM-labeled fusion-specific taqman assay (Thermo Fisher or IDT). qPCR assays were run on the QuantStudio 6 (Thermo Fisher, Waltham, MA) and Ct values were determined with the QuantStudio Real Time Software v1.3.

#### **RNA-Seq Data**

RNA-seq data was generated using the TruSeq RNA Exome kit (formerly RNA Access, Illumina, San Diego, CA) using 15 ng of total RNA as previously described (25). Libraries were sequenced on the NextSeq 500 platform (Illumina, San Diego, CA).

#### **RNA-Seq Data Analysis**

Fastq files were aligned to hg19 using STAR aligner, version 2.4.1b (32). Fusions were called using STAR-fusion version 0.5.4 (33). Variants were called using GATK version 3.3 (34), following the best practices for variant calling on RNA-seq (35, 36).

#### **Fusion Nomenclature**

All fusion partners are described in their 5'/3' order. In the commercial Afirma XA report, all fusion partners are reported alphabetically except CCDC6/RET and NCOA4/RET, which are reported using their colloquial names of RET/PTC1 and RET/PTC3, respectively.

#### Cohorts

943 blinded FNA samples with sufficient DNA were utilized from the following sources: Afirma GSC Algorithm Training set (n = 32, 217, 75, 35, 36 from Bethesda categories II-VI, respectively) (25), the prospectively collected, and multicenter Afirma GSC clinical validation cohort (n = 152 Bethesda III/IV, n = 17 Bethesda II and n = 29 Bethesda V/VI) (25), samples with paired castPCR *BRAF* V600E truth Bethesda V (n = 52) and Bethesda VI (n = 51) (37), and Bethesda III/IV samples from Afirma GEC (n = 247).

An additional 695 blinded FNAs from the Veracyte CLIA laboratory with sufficient RNA were deidentified and examined for fusions, during which time the total rate of assay failure among all samples received was 3.85%. This blinded and consecutive cohort was chosen without bias to represent an XA intended use cohort (n = 634 Bethesda III/IV GSC suspicious and n = 61 Bethesda V/VI).

#### **Statistics**

To evaluate the agreement when comparing the whole transcriptome RNA-seq to targeted AmpliSeq panels, positive percent agreement (PPA), negative percent agreement (NPA) and confirmation were calculated following the FDA guideline<sup>10</sup>. Statistical analyses were performed using R statistical software version 3.2.3<sup>11</sup>. All confidence intervals are 2-sided 95% CIs and were computed using the exact binomial test. The chi-square test of independence was performed to examine if there is a relationship between two categorial variables.

#### RESULTS

# A High Proportion of Variants Observed in DNA Are Expressed in RNA

To determine how many DNA variants could be detected in expressed RNA, RNA and DNA were first extracted from the same biological sample for direct comparison. Afirma GSC data was generated from the RNA, which utilizes whole-transcriptome RNA-seq data as an input to the machine learning algorithms. DNA was analyzed with a custom, targeted AmpliSeq panel that covers 761 variants that have been described in thyroid samples. Nine hundred forty-three FNA samples were analyzed with the custom DNA AmpliSeq panel. Four hundred forty-two samples were used for parameter tuning to ensure accurate detection of variants, and after the analysis pipeline was locked, the remaining 501 samples were used to evaluate the performance of calling a variant from RNA-seq data relative to DNA-based detection. Using a DNA variant allele frequency (VAF) cutoff of 5%, 181 DNA variants were observed, and the same variants were observed in 134 RNA-seq samples (74%; **Table 2**). Using a 20% VAF threshold for a positive result in DNA, positive percent agreement (PPA) increases to 88%.

#### Low Levels of Variant Allele Expression May Account for Some Differences in Variant Detection Between RNA and DNA Methods

To further investigate the role of transcription in the detection of expressed variants, a targeted variant panel using RNA as the template was employed rather than DNA. From the 943 FNA samples, a representative set of 102 FNAs that were variant positive by DNA AmpliSeq were tested for expressed variants with an RNA AmpliSeq variant panel. The PPA of RNA-seq whole-transcriptome variants rose from 76.5% vs. DNA in this subset to 88.9% vs. RNA (Table 2). Next, RNA Variant AmpliSeq data was examined for 17 samples that had a DNA variant identified, but no variant identified in the whole-transcriptome RNA-seq. Six of these 17 samples had dramatically different VAFs when comparing the DNA and RNA (Figure 1), with VAFs observed in DNA >10% while RNA-based VAFs were <5%. These six samples had a DNA:RNA VAF ratio ranging from 5.7 to 38.1 (including 3 samples with a DNA:RNA VAF ratio >10). In these six samples, the wild type allele is predominantly expressed and that biological difference accounts for the lack of variant detection in the whole-transcriptome RNA-seq data. In the most striking sample, we observed a VAF of 32% in the DNA and <1% in RNA.

#### There Is High Agreement Between Fusion Detection by Whole-Transcriptome RNA-Seq and Targeted Fusion Seq

To determine the PPA of the fusion-calling capabilities of wholetranscriptome RNA-seq data, an in-house custom RNA fusion AmpliSeq assay was developed. A new series of 695 consecutive and blinded FNAs from the Veracyte CLIA stream that were either Afirma GSC Suspicious or Bethesda V/VI were tested for fusions with the RNA fusion AmpliSeq assay. Any sample that generated a discordant result between whole-transcriptome RNA-seq and the RNA fusion AmpliSeq assay were resolved with qPCR. A total of 61 fusions were observed in this series, and the fusions observed in greater than one FNA were: PAX8/PPARG (n = 16), ETV6/NTRK3 (n = 13), RET/PTC1 (n = 6), STRN/ALK (n = 6), RET/PTC3 (n = 3), AGK/BRAF (n = 3), SND1/BRAF (n = 2), and *RBPMS/NTRK3* (n = 2) (**Supplementary Table 1**). This analysis revealed an 82.1% PPA between the RNA fusion AmpliSeq assay and whole-transcriptome RNA-seq data and demonstrated a 100% confirmation rate, as all 50 fusions

<sup>&</sup>lt;sup>10</sup>https://www.fda.gov/RegulatoryInformation/Guidances/ucm071148.htm section 7.2

<sup>11</sup> https://www.r-project.org

TABLE 2	Variant and Eusion	Performance in whol	e transcriptome RNA-se	a compared to targeter	d AmpliSeg panels
				g compared to targeter	

Genomic alteration	Samples	PPA	NPA	Confirmation	RNA-seq only	AmpliSeq only	Both detected
DNA Variants	501	74% [67–80]	100% [100–100]	98.5% [95–100]	2	47	134
RNA Variants	102	88.9% [80–95]	100% [100–100]	94.7% [87–99]	4	9	72
Fusions	695	82% [70–91]	100% [100–100]	100% [93–100]	0	11	50

Positive Percent Agreement (PPA) is percentage of variants or fusions that are positive by RNA-seq compared to AmpliSeq. Negative Percent Agreement (NPA) is the percentage of variants or fusions that are negative by RNA-seq compared to AmpliSeq. Confirmation is the percentage of variants or fusions that were called positive by RNA-seq that are also positive by AmpliSeq.





FIGURE 1 Variant Allele Frequencies (VAF) determined by targeted DNA and RNA AmpliSeq methods for 102 FNAs that were variant positive by DNA AmpliSeq. Samples with low RNA AmpliSeq coverage of the variant were excluded. Gray points were detected by all 3 methods, blue points were detected by DNA AmpliSeq and RNA AmpliSeq, but not detected by XA, red points were detected by DNA AmpliSeq but not RNA AmpliSeq or XA. The black dotted line is x=y.

identified by RNA-seq were also identified by the AmpliSeq assay (**Table 2**).

## Analytical Validation Shows the Variant and Fusion Calls Are Highly Reproducible Across Labs and Reagent Lots

To determine the reproducibility of the whole-transcriptome RNA-seq assay across laboratories and reagent lots, we examined the analytical validity of the assay. We compared the variant and fusion calls from RNA-seq data between two labs (R&D and CLIA) with the same lot of library prep reagents. Sixtynine variant positive samples were used, and the between lab accuracy was 89.9% (**Supplementary Figure 1A**). For fusions, 36 positive samples were used, and the between lab accuracy was 94.4% (**Supplementary Figure 1B**). Next, we examined the reproducibility of the assay within one plate and across different reagent lots. Nine variant positive samples and 6 fusion positive samples were plated in triplicate across 3 plates. Each plate was run with different reagent lots and different operators. These experiments investigated intra-plate and interplate reproducibility. For variants, the intra-plate reproducibility



was 88.9% and the inter-plate reproducibility was 86.4%. For fusions, the intra-plate reproducibility was 94.4% and the interplate reproducibility was 90.7%. These results passed prespecified acceptance criteria for these studies.

## Among Bethesda III/IV Nodules, 49% of Malignancies Harbor an RNA-Seq Detected Variant or Fusion

To understand the relationship between variants and fusions and histopathology diagnosis, we compared the variants/fusions observed in RNA-seq data from the primary test set of 190 Bethesda III/IV samples (Figure 2 and Table 3) that were collected in a prospective, multicenter, and blinded protocol for the clinical validation of the Afirma GEC (38), and subsequently utilized to clinically validate the Afirma GSC (25). The histopathological diagnosis of these nodules was assigned by an expert panel of thyroid histopathologists who were masked to all clinical, molecular, and cytological data. In 145 histologically benign nodules, 76% had no variant or fusion observed (76% specificity). Of the 24% with a variant, RAS variants were the most common, followed by TSHR. In the 45 histologically malignant samples, 51% had no variant or fusion observed (49% sensitivity). Overall, the RNA-seq PPV and NPV were 38% and 83%, respectively, in this cohort with a 24% cancer prevalence. The 49% of samples with a positive variant or fusion mostly harbored RAS variants, BRAF variants, or fusions. Of the variants observed more than once, only BRAF V600E was confined to malignant nodules. Conversely, 10 of 11 TSHR variants and 4 of 5 SPOP variants occurred in histologically benign nodules. The only gene fusion observed more than once was *PAX8/PPARG* and both occurrences were in histologically benign nodules. *MKRN1/BRAF* and *ETV6/NTRK3* fusions were each identified once and were in malignant nodules (both PTC). *PAX8/GLIS3* was identified in the one hyalinizing trabecular adenoma, consistent with a recent report (39).

We next examined the combination of Afirma GSC classifier prediction and variant/fusion results (Figure 2 and Table 3). In the 99 histologically benign nodules with GSC benign results (true negatives), 15% were variant positive and none contained a fusion. In the 46 benign nodules with GSC suspicious results (false positives), 39% harbored a variant and 6.5% contained a fusion (2 PAX8/PPARG and 1 PAX8/GLIS3). In the 41 histologically malignant nodules with GSC suspicious results (true positives), 46% had a variant, and 5% had a fusion (MKRN1/BRAF and ETV6/NTRK3). Finally, in the 4 GSC benign false negative nodules (2 PTC, 1 FVPTC, 1 HCC), only the HCC contained a variant (TSHR). Taken together, in the 190 thyroid nodules with definitive histology, malignant nodules were twice as likely to carry a variant or fusion relative to benign nodules (49 vs. 24%, p = 0.003 [y 2]). The most common genomic alteration identified among GSC suspicious nodules was a variant in a RAS family gene (present in 32% and conveying a PPV of 46%). Taken together, GSC suspicious nodules carried a variant or fusion in 48% compared with 16% of GSC benign nodules (p < 0.0001). The most common variants among GSC benign nodules were TSHR, and no variants or fusions known to be highly predictive of malignancy were identified among nodules with GSC benign results. In clinical use, variants and fusions are not reported among GSC benign nodules.

Among GSC suspicious nodules with a cancer prevalence of 47%, the sensitivity, specificity, and NPV of XA were 51,


			XA variant and fusions (n)		
Histopathological subtype	Nodules, No. (%)	GSC B No./GSC S No.	GSC benign	GSC suspicious	
BENIGN					
Total	145 (100)	99/46			
BFN	49 (33.8)	38/11	EIF1AX:p.G8R(1)	HRAS:p.Q61R(2)	
			TSHR:p.I486M(1)	NRAS:p.Q61K(1)	
			TSHR:p.L512R(1)	NRAS:p.Q61R(1)	
			TSHR:p.M453T(1)		
HN	5 (3.4)	5/0	TSHR:p.I568T(1)	NA	
FA	54 (37.2)	37/17	GNAS:p.Q870H(1)	HRAS:p.Q61R(1)	
			SPOP:p.P94R(3)	NRAS:p.Q61R(3)	
			TSHR:p.D633Y(1)	PAX8/PPARG(1)	
			TSHR:p.L629F(1)		
-T-UMP	9 (6.2)	4/5	NA	HRAS:p.Q61R(1)	
				NRAS:p.Q61R(1)	
				SPOP:p.P94R(1)	
NDT-UMP	8 (5.5)	4/4	NA	HRAS:p.Q61R(2)	
				NRAS:p.Q61R(1)	
				PAX8/PPARG(1)	
HCA	17 (11.7)	10/7	EIF1AX:p.G9D(1)	NRAS:p.Q61K(2)	
			PTEN:p.G129R(1)	TSHR:p.M453T(1)	
			TSHR:p.L629F, EZH1:p.Y642F(1)	TSHR:p.S281I(1)	
			TSHR:p.S425I(1)		
CLT	2 (1.4)	1/1	NA	NA	
HTA	1 (0.7)	0/1	NA	PAX8/GLIS3(1)	
MALIGNANT					
Total	45 (100)	4/41			
PTC	15 (33.3)	2/13	NA	BRAF:p.V600E(3)	
				NRAS:p.Q61R(1)	
				SPOP:p.P94R(1)	
				MKRN1/BRAF (1)	
				ETV6/NTRK3(1)	
PTC-TCV	1 (2.2)	0/1	NA	NA	
FV-PTC	11 (24.4)	1/10	NA	HRAS:p.Q61R(1)	
				KRAS:p.Q61R, EIF1AX:p.A113_splice(1)	
				NRAS:p.Q61K(1)	
				NRAS:p.Q61R(4)	
HCC-c	9 (20)	1/8	TSHR:p.I568T(1)	EIF1AX:p.A113_splice(1)	
				NRAS:p.Q61R(1)	
FC	7 (15.6)	0/7	NA	BRAF:p.K601E(1)	
	· · ·			HRAS:p.G13R(1)	
				NRAS:p.Q61R(1)	
PDTC	1 (2.2)	0/1	NA	NRAS:p.Q61K(1)	

TABLE 3 | Variants and fusions relative to their histopathologic and Afirma GSC outcomes (25). GSC Benign (B), GSC Suspicious (S).

54, and 56%, respectively. The overall PPV of an Afirma GSC suspicious nodule was 47%, regardless of variant/fusion status. The PPV was 50% among GSC suspicious nodules when a variant or fusions was identified, compared with 44% among GSC suspicious nodules when no variant or fusion was identified (p = 0.77 [ $\chi$ 2]).

#### Variants and Fusions Potentially Amenable to Targeted Therapy Were Identified Across Multiple Bethesda Categories

We determined the most common variants observed in thyroid FNAs from expressed variant data (Figure 3 and Supplementary Table 2). Five genes showed variants in >1% of



FNAs: *BRAF*, *NRAS*, *HRAS*, *TSHR*, and *SPOP*. For individual variants, *BRAF* V600E was most common, followed by *NRAS* Q61R, *HRAS* Q61R, *NRAS* Q61K, *TSHR* M453T, and *SPOP* P94R (**Supplementary Table 2**). *NRAS* and *HRAS* variants were primarily observed in Bethesda III/IV FNAs that were GSC Suspicious, while *BRAF* V600E was predominantly present in Bethesda V/VI FNAs. *TSHR* and *SPOP* variants were most frequently observed in GSC Benign FNAs (**Figure 3**), which is consistent with other observations (11, 40).

Fusions were also identified using whole-transcriptome RNA-seq. Fifty-two (5.5%) of the 943 samples harbored a fusion. The most common fusions were *CCDC6/RET* (aka *RET/PTC1*), *PAX8/PPARG*, and *ETV6/NTRK3* (Figure 3 and Supplementary Table 3). *RET/PTC1* was primarily observed in Bethesda V/VI FNAs, while *PAX8/PPARG* was observed in Bethesda III/IV FNAs. *ETV6/NTRK3* was observed across Bethesda III-VI. A subset of 31 fusion-containing samples with sufficient RNA were selected for testing by qPCR (Supplementary Table 4). 100% of the RNA-seq detected fusions were also detected (i.e., confirmed) by qPCR. *ALK* fusions were also observed, although they were rare. *EML4/ALK* was observed twice in Bethesda VI samples and *STRN/ALK* was observed once in Bethesda III.

Overall, variants or fusions potentially amenable to targeted therapy involving *AKT*, *ALK*, *BRAF*, *HRAS*, *KRAS*, *MET*, *NRAS*, *NTRK*, *PPARG*, *PTEN*, and *RET* were identified across multiple Bethesda categories. Similarly, 4 *RET* point mutations suggestive of MTC were identified among Bethesda III-V samples (**Supplementary Table 2**). All 4 were identified as positive by the Afirma GSC MTC classifier (25, 41). Further, these variants may be somatic or germline and raise the possibility of multiple endocrine neoplasia type 2.

## 1% of Bethesda III/IV Nodules Harbor a TERT Promoter Variant, All Were in Combination With a RAS Variant and Were GSC Suspicious, and Most Were Histologically Benign

The DNA AmpliSeq panel included the TERT promoter. TERT promoter variants were observed in 15 samples of the 943 examined, across Bethesda categories (Table 4), and all were called Afirma GSC suspicious: 1 Bethesda II (1.89%), 0 Bethesda III (0%), 7 Bethesda IV (3.48%), 2 Bethesda V (1.90%), and 6 Bethesda VI (5.61%). The Bethesda II sample was TERT C228T plus NRAS Q61K. All 7 Bethesda IV FNAs were TERT C228T plus RAS positive. Seven of eight Bethesda V and VI FNAs were TERT C228T plus BRAF V600E, with the remaining Bethesda VI sample TERT C228T positive in isolation. Fourteen of fifteen TERT C228T positive FNAs had paired histopathology (Table 4). The Bethesda II sample was an FVPTC. In Bethesda IV, 5 FNAs were TERT C228T plus NRAS positive (4 Q61R and 1 Q61K) and 4/5 were histologically benign. Additionally, there were 2 FNAs that were positive for TERT C228T in combination with KRAS Q61R or HRAS Q61R, which were histologically malignant and benign, respectively. Finally, five Bethesda V and VI FNAs were TERT C228T in combination with BRAF V600E, and all 5 were histologically malignant. Three were PTC and two were PTC-TCV. The Bethesda VI sample with TERT C228T in isolation was an FC-v.

TABLE 4	FERT Promoter variants observed in this st	tudy.
		uuy.

Bethesda	Histopathology	SubType	Afirma GSC	Variant	TERT Promoter Variant
Bethesda II	М	FVPTC	Suspicious	NRAS:p.Q61K	C228T
Bethesda IV	В	HCA	Suspicious	HRAS:p.Q61R	C228T
	В	HCA	Suspicious	NRAS:p.Q61K	C228T
	В	NHP	Suspicious	NRAS:p.Q61R	C228T
	В	FA	Suspicious	NRAS:p.Q61R	C228T
	В	WDT-UMP	Suspicious	NRAS:p.Q61R	C228T
	М	PTC	Suspicious	KRAS:p.Q61R	C228T
	М	FC-v	Suspicious	NRAS:p.Q61R	C228T
Bethesda V	М	PTC	Suspicious	BRAF:p.V600E	C228T
	М	PTC	Suspicious	BRAF:p.V600E	C228T
Bethesda VI	М	PTC	Suspicious	BRAF:p.V600E	C228T
	М	PTC-TCV	Suspicious	BRAF:p.V600E	C228T
	Μ	PTC-TCV	Suspicious	BRAF:p.V600E	C228T
	Unknown	Unknown	Suspicious	BRAF:p.V600E	C228T
	Μ	FC-v	Suspicious	None	C228T

Histopathology is shown as Malignant (M) or Benign (B). See Table 1 for a list of subtype abbreviations. Unknown subtype has no histopathology information available.

#### DISCUSSION

This study demonstrates the analytical and clinical validation in thyroid nodule evaluation of the Afirma Xpression Atlas (XA), which detects gene variants and fusions from a curated panel of 511 genes via dedicated FNA samples using whole-transcriptome RNA-sequencing. More than 96% of consecutive real-world samples received has sufficient quantity and quality to receive an XA result. The data show that agreement is high between fusion detection by XA and an alternative fusion detection method, and that a high proportion of variants observed in DNA are also detected in the expressed RNA. As expected, agreement is even higher when comparing variant detection between two different RNAbased methods. Intriguingly, we found the wild type allele to be preferentially expressed compared with the variant allele in some samples (Figure 1), explaining some of the differences observed. Analytical studies show high reproducibility within plates and across reagent lots. Among Bethesda III/IV cytologically indeterminate nodules that were malignant, about half harbored a variant or fusion that was detected by XA. Detection of variants and fusions progressively increased along the Bethesda II to VI spectrum, and genomic findings potentially amenable to targeted therapeutics were identified across the Bethesda spectrum.

Afirma XA is not a cancer rule-out test. Among Bethesda III/IV nodules deemed suspicious by Afirma GSC and among Bethesda V/VI nodules, the impact of Afirma XA on nodule management extends beyond informing the risk of cancer when XA is negative or when XA is positive for a specific genomic alteration. Recent studies have begun to associate selected variant and fusions with *BRAF* V600E-like vs. *RAS*-like (or non-BRAF-non-RAS) pathway signaling, iodine metabolism, neoplasm histology, risk of lymph node metastasis, risk of recurrence, and risk of mortality (3, 9). How much these

prognostic associations will remain significant independent predictors when traditional predictors are considered is presently unknown. For instance, TERT promoter mutations predict disease-free survival and disease-specific survival, but this effect is diminished or eliminated when the variant occurs in the absence of a RAS or BRAF variant, or among low and intermediate ATA risk patients and stage I-II TNM patients (42). Thus, for thyroid nodules seemingly confined to the thyroid gland and harboring Bethesda III-VI cytology, future studies may investigate how the preoperative identification of a presumed driver mutation may inform the pre-operative evaluation and surgical plan. For example, prospective randomized trials could investigate the roles of active surveillance or hemithyroidectomy for variants and fusions associated with a ~50:50 chance of cancer, and/or those associated with less aggressive carcinomas. This may include nodules where a variant or fusion is not detected as some data suggests that such cancers may have less aggressive features when considering extrathyroidal extension, lymph node metastases, risk of recurrence, and risk of mortality (3, 9). Similarly, prospective randomized trials could investigate the role of total thyroidectomy vs. hemithyroidectomy (or whether or not adjuvant radioactive iodine ablation or TSH suppression would be used) for cancers with variants thought to predict more aggressive tumor behavior that are nevertheless clinically confined to the thyroid gland.

For patients presenting with locally advanced thyroid cancer, genomic findings on XA may suggest the possibility of neoadjuvant therapy that may improve the outcome of a subsequent surgical resection (16, 43, 44).

For patients with thyroid carcinoma in the neck or in distant sites that is refractory to radioactive iodine and may warrant systemic therapy, use of XA from FNAs of known thyroid cancer deposits may inform treatment selection, although confirmation testing by an approved companion diagnostic test may still be required for patients to access certain pharmaceuticals. Additionally, repeat testing from a site of disease progression during active treatment may provide additional genomic insights to potentially guide therapy (26).

Limitations of measuring variants in expressed RNA include that some variants and fusions identified by an alternative method were not identified by XA. The reason for these differences is not known, nor is it known which test should ideally be considered "correct." While imperfect test sensitivity is one possibility, it is also possible that some DNA variants may not be expressed due to gene silencing, or very low expression levels. Such phenomenon may explain some discrepancies among BRAF V600E variants detected by qPCR that are negative by immunohistochemistry (45). By employing a third variant detection methodology that used targeted sequencing of RNA templates, we have shown that some samples do have very low expression of the gene variants identified in DNA. The biological significance of such variants is unknown. The efficacy of targeted treatment aimed at non-expressed or poorly expressed genomic alterations may be diminished. Conversely, the vast majority of genomic abnormalities identified by XA were confirmed by the alternative method (Table 2). An additional limitation of RNA sequencing is that variants in non-coding regions, such as TERT promotor variants, are not detected by this method. However, our data demonstrate that these variants are uncommon among cytologically indeterminate nodules (<1%) and in the vast majority of cases, found in tandem with a RAS variant. Moreover, we show that TERT promoter variants in combination with RAS variants can occur in benign lesions in Bethesda IV FNAs. While current opinion is that nodules with a RAS variant (with or without TERT promoter variant) should be surgically removed given their potential malignant or pre-malignant status, it is unclear if cancers harboring a TERT promoter variant plus a RAS or BRAF variant should be treated differently based on this genomic information independent from traditional prognostic factors for risks of recurrence and death, especially among lower-risk patients (42). A DNA based detection method, or development of an RNA expression-based classifier, could be added to XA in the future should reporting of such nonexpressed variants be desired.

In clinical practice, XA testing is offered for Bethesda III/IV Afirma GSC suspicious nodules, FNA Bethesda V/VI nodules, and for known thyroid cancer metastases. The Afirma Xpression Atlas became commercially available to most of the United States in 2018. Samples for Afirma GSC and/or XA are collected with two dedicated FNA needle passes that must be expressed into the provided FNAprotect tube, properly stored locally to avoid exposure to heat, and are shipped in the provided container with frozen foam bricks to ensure receipt of high-quality RNA material. These validated collection and shipping procedures are identical to those used to formerly collect and ship Afirma GEC samples (46). Detailed sample collection, packing, and shipping instructions are available on-line<sup>12</sup>.

In summary, we have demonstrated clinical and analytical validation of the Afirma XA, which reports variants and fusions from a panel of 511 genes that have been associated with

thyroid cancer. This added clinical information is intended to supplement clinical decision-making among patients with Bethesda III-VI nodules. Clinicians are to be reminded that most patients with thyroid cancer have an excellent prognosis, and the greatest impact of this added genomic information may be to facilitate treatment that is less aggressive, rather than more aggressive (47). The information obtained from variants and fusions assessment may offer new precision medicine insights from diagnostic FNA samples and the opportunity to advance individualized patient care.

## DATA AVAILABILITY

Restrictions apply to the datasets: The datasets for this manuscript are not publicly available because the dataset and the research methodologies are proprietary. Requests to access the datasets should be directed to GK, Ph.D.; Giulia@Veracyte.com.

#### ETHICS STATEMENT

This research was approved by the Copernicus Group Independent Review Board (Cary, North Carolina). A waiver of written informed consent was granted regarding de-identified biological materials from the CLIA laboratory. IRB approval and written Informed consent in accordance with the Declaration of Helsinki was provided by all patients whose samples were previously used for training and validation of the Afirma GSC as previously described (25, 38).

## **AUTHOR CONTRIBUTIONS**

GK, JB, JH, PW, SK, and YH conceived and designed the study. JB, LW, PW, RK, SK, TA, and YH analyzed and interpreted the data. EC, GK, JB, LW, MC, MS, RK, SW, and TA drafted and critically revised the work for important intellectual content. SK and YH performed the bioinformatics and statistical analyses. All authors contributed to manuscript revision, read and approved the submitted version.

#### **FUNDING**

This research was funded by Veracyte, Inc.

## ACKNOWLEDGMENTS

We would like to thank Ed Tom, Duncan Whitney, Zhanzhi Hu, Mei Wong, Grazyna Fedorowicz, Manqiu Cao, Huimin Jiang, Jessica Anderson, Tami Tu, Maggie Dorosz, Jennifer Huiras, and Jing Lu for their contributions to laboratory studies and data analysis.

#### SUPPLEMENTARY MATERIAL

The Supplementary Material for this article can be found online at: https://www.frontiersin.org/articles/10.3389/fendo. 2019.00612/full#supplementary-material

<sup>&</sup>lt;sup>12</sup>www.afirma.com/physicians/practice-resources

## REFERENCES

- Stratton MR, Campbell PJ, Futreal PA. The cancer genome. Nature. (2009) 458:719–24. doi: 10.1038/nature07943
- Agrawal N, Jiao Y, Sausen M, Leary R, Bettegowda C, Roberts NJ, et al. Exomic sequencing of medullary thyroid cancer reveals dominant and mutually exclusive oncogenic mutations in RET and RAS. J Clin Endocrinol Metab. (2013) 98:E364–369. doi: 10.1210/jc.2012-2703
- Cancer Genome Atlas Research Network. Integrated genomic characterization of papillary thyroid carcinoma. *Cell.* (2014) 159:676–90. doi: 10.1016/j.cell.2014.09.050
- Costa V, Esposito R, Ziviello C, Sepe R, Bim LV, Cacciola NA, et al. New somatic mutations and WNK1-B4GALNT3 gene fusion in papillary thyroid carcinoma. *Oncotarget*. (2015) 6:11242–51. doi: 10.18632/oncotarget.3593
- Jung SH, Kim MS, Jung CK, Park HC, Kim SY, Liu J, et al. Mutational burdens and evolutionary ages of thyroid follicular adenoma are comparable to those of follicular carcinoma. *Oncotarget*. (2016) 7:69638–48. doi: 10.18632/oncotarget.11922
- Pagan M, Kloos RT, Lin CF, Travers KJ, Matsuzaki H, Tom EY, et al. The diagnostic application of RNA sequencing in patients with thyroid cancer: an analysis of 851 variants and 133 fusions in 524 genes. *BMC Bioinform*. (2016) 17(Suppl. 1):6. doi: 10.1186/s12859-015-0849-9
- Pan W, Zhou L, Ge M, Zhang B, Yang X, Xiong X, et al. Whole exome sequencing identifies lncRNA GAS8-AS1 and LPAR4 as novel papillary thyroid carcinoma driver alternations. *Hum Mol Genet.* (2016) 25:1875–84. doi: 10.1093/hmg/ddw056
- Siraj AK, Masoodi T, Bu R, Beg S, Al-Sobhi SS, Al-Dayel F, et al. Genomic profiling of thyroid cancer reveals a role for thyroglobulin in metastasis. *Am J Hum Genet.* (2016) 98:1170–80. doi: 10.1016/j.ajhg.2016.04.014
- Yoo SK, Lee S, Kim SJ, Jee HG, Kim BA, Cho H, et al. Comprehensive analysis of the transcriptional and mutational landscape of follicular and papillary thyroid cancers. *PLoS Genet.* (2016) 12:e1006239. doi: 10.1371/journal.pgen.1006239
- Lu Z, Zhang Y, Feng D, Sheng J, Yang W, Liu B. Targeted next generation sequencing identifies somatic mutations and gene fusions in papillary thyroid carcinoma. *Oncotarget*. (2017) 8:45784–92. doi: 10.18632/oncotarget.17412
- Ye L, Zhou X, Huang F, Wang W, Qi Y, Xu H, et al. The genetic landscape of benign thyroid nodules revealed by whole exome and transcriptome sequencing. *Nat Commun.* (2017) 8:15533. doi: 10.1038/ncomms15533
- Ganly I, Makarov V, Deraje S, Dong Y, Reznik E, Seshan V, et al. Integrated genomic analysis of hurthle cell cancer reveals oncogenic drivers, recurrent mitochondrial mutations, and unique chromosomal landscapes. *Cancer Cell*. (2018) 34:256–70 e255. doi: 10.1016/j.ccell.2018.07.002
- Gopal RK, Kubler K, Calvo SE, Polak P, Livitz D, Rosebrock D, et al. Widespread chromosomal losses and mitochondrial DNA alterations as genetic drivers in hurthle cell carcinoma. *Cancer Cell.* (2018) 34:242–55 e245. doi: 10.1016/j.ccell.2018.06.013
- Pozdeyev N, Gay LM, Sokol ES, Hartmaier R, Deaver KE, Davis S, et al. Genetic analysis of 779 advanced differentiated and anaplastic thyroid cancers. *Clin Cancer Res.* (2018) 24:3059–68. doi: 10.1158/1078-0432.CCR-18-0373
- Drilon A, Siena S, Ou SI, Patel M, Ahn MJ, Lee J, et al. Safety and antitumor activity of the multitargeted pan-TRK, ROS1, and ALK inhibitor entrectinib: combined results from two phase I trials (ALKA-372-001 and STARTRK-1). *Cancer Discov.* (2017) 7:400–9. doi: 10.1158/2159-8290.CD-16-1237
- Rao SN, Cabanillas ME. Navigating systemic therapy in advanced thyroid carcinoma: from standard of care to personalized therapy and beyond. J Endocr Soc. (2018) 2:1109–30. doi: 10.1210/js.2018-00180
- Subbiah V, Kreitman RJ, Wainberg ZA, Cho JY, Schellens JHM, Soria JC, et al. Dabrafenib and trametinib treatment in patients with locally advanced or metastatic BRAF V600-mutant anaplastic thyroid cancer. *J Clin Oncol.* (2018) 36:7–13. doi: 10.1200/JCO.2017.73.6785
- Untch BR, Dos Anjos V, Garcia-Rendueles MER, Knauf JA, Krishnamoorthy GP, Saqcena M, et al. Tipifarnib inhibits HRASdriven dedifferentiated thyroid cancers. *Cancer Res.* (2018) 78:4642–57. doi: 10.1158/0008-5472.CAN-17-1925
- Tian T, Li X, Zhang J. mTOR signaling in cancer and mTOR inhibitors in solid tumor targeting therapy. *Int J Mol Sci.* (2019) 20:E755. doi: 10.3390/ijms20030755

- Drilon A, Laetsch TW, Kummar S, DuBois SG, Lassen UN, Demetri GD, et al. Efficacy of larotrectinib in TRK fusion-positive cancers in adults and children. N Engl J Med. (2018) 378:731–9. doi: 10.1056/NEJMoa17 14448
- Thorpe LM, Yuzugullu H, Zhao JJ. PI3K in cancer: divergent roles of isoforms, modes of activation and therapeutic targeting. *Nat Rev Cancer*. (2015) 15:7–24. doi: 10.1038/nrc3860
- Hu M, Taylor M, Wirth LJ, Zhu VW, Doebele RC, Lee D, et al. Clinical activity of selective RET inhibitor, BLU-667, in advanced RET-altered thyroid cancers: updated results from the phage 1 ARROW study (short call oral 5). *Thyroid.* (2018) 28, A-170.
- Wirth LJ, Cabanillas ME, Sherman E, Solomon B, Leboulleux S, Robinson B, et al. Clinical activity of LOXO-292, a highly selective RET inhibitor, in patients with RET-altered thyroid cancers (short call oral 6). *Thyroid*. (2018) 28, A-171.
- 24. Marcus L, Lemery SJ, Keegan P, Pazdur R. FDA approval summary: pembrolizumab for the treatment of microsatellite instability-high solid tumors. *Clin Cancer Res.* (2019) 25:3753–8. doi: 10.1158/1078-0432.CCR-18-4070
- Patel KN, Angell TE, Babiarz J, Barth NM, Blevins T, Duh QY, et al. Performance of a genomic sequencing classifier for the preoperative diagnosis of cytologically indeterminate thyroid nodules. *JAMA Surg.* (2018) 153:817– 24. doi: 10.1001/jamasurg.2018.1153
- 26. Ali SZ, Siperstein A, Sadow PM, Golding AC, Kennedy GC, Kloos RT, et al. Extending expressed RNA genomics from surgical decision making for cytologically indeterminate thyroid nodules to targeting therapies for metastatic thyroid cancer. *Cancer Cytopathol.* (2019). doi: 10.1002/cncy.22132. [Epub ahead of print].
- 27. Cibas ES, Ali SZ. The 2017 bethesda system for reporting Thyroid Cytopathology. *Thyroid*. (2017) 27:1341-6. doi: 10.1089/thy.2017.0500
- Parkinson D, Aziz S, Bentley R, Johnson SJ. Thyroid cytology-histology correlation using the RCPath terminology for thyroid cytology reporting. J Clin Pathol. (2017) 70:648–55. doi: 10.1136/jclinpath-2016-204022
- Angell TE, Heller HT, Cibas ES, Barletta JA, Kim MI, Krane JF, et al. Independent comparison of the afirma genomic sequencing classifier and gene expression classifier for cytologically indeterminate Thyroid Nodules. *Thyroid*. (2019) 29:650–6. doi: 10.1089/thy.2018.0726
- Endo M, Nabhan F, Porter K, Roll K, Shirley L, Azaryan I, et al. Afirma gene sequencing classifier compared to gene expression classifier in indeterminate Thyroid Nodules. *Thyroid*. (2019) 29:1115–24. doi: 10.1089/thy.2018.0733
- Harrell RM, Eyerly-Webb SA, Golding AC, Edwards CM, Bimston DN. Statistical comparison of afirma gsc and afirma gec outcomes in a community endocrine surgical practice: early findings. *Endocr Pract.* (2019) 25:161–4. doi: 10.4158/EP-2018-0395
- Dobin A, Davis CA, Schlesinger F, Drenkow J, Zaleski C, Jha S, et al. STAR: ultrafast universal RNA-seq aligner. *Bioinformatics*. (2013) 29:15–21. doi: 10.1093/bioinformatics/bts635
- Haas B, Dobin A, Stransky N, Li B, Yang X, Tickle T, et al. STAR-Fusion: fast and accurate fusion transcript detection from RNA-Seq. *bioRxiv*. (2017). doi: 10.1101/120295
- McKenna A, Hanna M, Banks E, Sivachenko A, Cibulskis K, Kernytsky A, et al. The Genome Analysis Toolkit: a MapReduce framework for analyzing next-generation DNA sequencing data. *Genome Res.* (2010) 20:1297–303. doi: 10.1101/gr.107524.110
- 35. DePristo MA, Banks E, Poplin R, Garimella KV, Maguire JR, Hartl C, et al. A framework for variation discovery and genotyping using next-generation DNA sequencing data. *Nat Genet*. (2011) 43:491–8. doi: 10.1038/ng.806
- Van der Auwera GA, Carneiro MO, Hartl C, Poplin R, Del Angel G, Levy-Moonshine A, et al. From FastQ data to high confidence variant calls: the Genome Analysis Toolkit best practices pipeline. *Curr Protoc Bioinformatics*. (2013) 43, 11–33. doi: 10.1002/0471250953.bi1110s43
- Diggans J, Kim SY, Hu Z, Pankratz D, Wong M, Reynolds J, et al. Machine learning from concept to clinic: reliable detection of BRAF V600E DNA mutations in thyroid nodules using high-dimensional RNA expression data. *Pac Symp Biocomput*. (2015) 371–82.
- Alexander EK, Kennedy GC, Baloch ZW, Cibas ES, Chudova D, Diggans J, et al. Preoperative diagnosis of benign thyroid nodules with indeterminate cytology. N Engl J Med. (2012) 367:705–15. doi: 10.1056/NEJMoa1203208

- Nikiforova MN, Nikitski AV, Panebianco F, Kaya C, Yip L, Williams M, et al. GLIS rearrangement is a genomic hallmark of hyalinizing trabecular tumor of the *Thyroid* Gland. *Thyroid*. (2019) 29:161–73. doi: 10.1089/thy. 2018.0791
- Guan H, Matonis D, Toraldo G, Lee SL. Clinical significance of thyroidstimulating hormone receptor gene mutations and/or sodium-iodine symporter gene overexpression in indeterminate thyroid fine needle biopsies. *Front Endocrinol.* (2018) 9:566. doi: 10.3389/fendo.2018.00566
- 41. Randolph G, Angell TE, Babiarz J, Barth N, Blevins T, Duh Q, et al. Clinical validation of the afirma genomic sequencing classifier for medullary thyroid cancer (Clinical Oral Abstract 29). *Thyroid*. (2017) 27:A105.
- 42. Song YS, Lim JA, Choi H, Won JK, Moon JH, Cho SW, et al. Prognostic effects of TERT promoter mutations are enhanced by coexistence with BRAF or RAS mutations and strengthen the risk prediction by the ATA or TNM staging system in differentiated thyroid cancer patients. *Cancer*. (2016) 122:1370–9. doi: 10.1002/cncr.29934
- Cabanillas ME, Busaidy NL, Zafereo M, Waguespack SG, Hu MI, Hofmann MC, et al. Neoadjuvant vemurafenib in patients with locally advanced papillary thyroid cancer (abstract). *Eur Thyroid J.* 6(supp 1):38. doi: 10.1159/000477987
- 44. Cabanillas ME, Ferrarotto R, Garden AS, Ahmed S, Busaidy NL, Dadu R, et al. Neoadjuvant BRAF- and immune-directed therapy for anaplastic *Thyroid* carcinoma. *Thyroid*. (2018) 28:945–51. doi: 10.1089/thy.20 18.0060
- 45. Szymonek M, Kowalik A, Kopczynski J, Gasior-Perczak D, Palyga I, Walczyk A, et al. Immunohistochemistry cannot replace DNA analysis for evaluation

of BRAF V600E mutations in papillary thyroid carcinoma. *Oncotarget*. (2017) 8:74897–909. doi: 10.18632/oncotarget.20451

- Walsh PS, Wilde JI, Tom EY, Reynolds JD, Chen DC, Chudova DI, et al. Analytical performance verification of a molecular diagnostic for cytologyindeterminate thyroid nodules. J Clin Endocrinol Metab. (2012) 97:E2297– 2306. doi: 10.1210/jc.2012-1923
- Valderrabano P, Khazai L, Thompson ZJ, Leon ME, Otto KJ, Hallanger-Johnson JE, et al. Impact of oncogene panel results on surgical management of cytologically indeterminate thyroid nodules. *Head Neck*. (2018) 40:1812–23. doi: 10.1002/hed.25165

**Conflict of Interest Statement:** GK, JB, JH, PW, RK, SK and YH are Veracyte Inc. employees and equity owners.

The remaining authors declare that the research was conducted in the absence of any commercial or financial relationships that could be construed as a potential conflict of interest.

Copyright © 2019 Angell, Wirth, Cabanillas, Shindo, Cibas, Babiarz, Hao, Kim, Walsh, Huang, Kloos, Kennedy and Waguespack. This is an open-access article distributed under the terms of the Creative Commons Attribution License (CC BY). The use, distribution or reproduction in other forums is permitted, provided the original author(s) and the copyright owner(s) are credited and that the original publication in this journal is cited, in accordance with accepted academic practice. No use, distribution or reproduction is permitted which does not comply with these terms.