

**Figure 1.**The Laerdal Pocket Infant Resuscitator Mask (participant's permission obtained)



**1a.** A close up of the Laerdal mask. The re-usable mask is fitted with a soft infant-sized silicone rubber face cup. The breathing tube is offset to allow for better observation of the infant during assisted ventilation. It has a one-way valve, preventing oral secretions from being aspirated by the rescuer. TBAs were trained to clean it after each use using soap and water.

**1b.** The small size of the mask and lack of a reservoir air bag facilitates its use by a solo rescuer, as demonstrated here.